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Analytical technique: pyrolysis-gas chromatography/mass spectrometry (py-GC/MS)

Sampling and analysis:

suitable for qualitative and semi-quantitative analysis of organic compounds with particular attention to polymers, in liquid and solid samples. Samples may be analyzed without pretreatment, with direct insertion of a small amount of sample (~ 0,1 mg) in a quartz tube, with or without derivatizing agents. The derivatizant may be added also in a preliminary step, followed by warming to obtain enhancement of hydrolysis reactions. The analysis is microinvasive and microdestructive.

Functioning principle:

The sample is loaded in a small quartz tube, retained with quartz wool and inserted in the pyrolyzer coil (probe). The probe is then inserted in a thermostated pyrolysis chamber, directly interfaced with GC/MS system. High temperature pyrolysis is started simultaneously with the chromatographic run, so chemical specie generated from pyrolysis are pushed by helium flux in the thermostated column, separated and relaved by mass spectrometer.

Analytical procedure for ancient ointment samples:

Two different procedures were followed for sample preparation. In the first one a small amount of sample (~ 0,1 mg) was loaded in the quartz tube without any pretreatment and 5/6 • l of TMAH aq. 25 % w/w were added before analysis. In the second one the sample (1-2 mg) was dissolved in TMAH aqueous solution (50-100 • l). The solution was submitted to warming in oil bath at 100 °C for 30 minutes, then 5-7 • l were inserted in the quartz tube and pyrolyzed.

Instrumentation and experimental conditions:

A Gaschromatograph model 6890N (Agilent Technologies, USA) coupled with mass spectrometer 5973 Mass Selective Detector (Agilent Technologies, USA) with quadrupole and electronic impact ionization (70eV) was used. The Pyrolyzer interfaced with the GC was a CDS Pyroprobe 1000 with CDS 1500 Valved Interface (CDS ANALYTICAL INC, USA). Pyrolysis and GC/MS interfaces were kept at 280 °C. Pyrolysis temperature was 600 °C.

GC column is an HP-5MS (Hewlett Packard, stationary phase 5% phenylmethylsiloxane, i.d. 0,25 mm, stationary phase thickness 0,25 • m, length 30 m). Separation conditions were:

- Oven program: 50 °C (held 2 min.), 10°C/min to 130 °C, 5 °C/min to 180 °C, 15°C/min to 300 °C (held 5 min.)
- Carrier gas : helium at 1,0 ml/min
- Injection in split mode at 1720 of total flux


Informations obtained:

Py-GC/MS analysis allows to obtain the molecular profile of organic compounds derived from pyrolysis of natural and synthetic polymeric materials. The kind of markers and their distribution make possible to recognize materials as natural and synthetic resins, polysaccharides, siccative oils and other polymers. It is also possible to obtain molecular profile of non-polymeric compounds, as oils and waxes, composed by mixtures of high molecular weight organic molecules (long chain fatty acids, alcohols and alkanes). More difficult is the proteinaceous material identification, because of the strong degradation of polypeptides at the high temperatures of analysis.

In addition to the characterization of substances achievable by identification of the specific markers formed during the pyrolysis, semi-quantitative compositional analysis may be obtained from peak areas evaluation.

Example of analysis:

Pyrogram of sample “*unguentum pro igne*” from Museo Aboca (San Sepolcro, Arezzo, Italy), n° 50015. Solid sampled charged in quartz tube with additions of 7 • l of TMAH aq. at 25 % w/w.

Saturated fatty acids methyl esters of C16 and C18 (palmitic and stearic acids) and unsaturated fatty acids methyl esters C18 (oleic acids) are indicative of the presence of plant oil, also confirmed by the peaks of C8 and C9 dicarboxylic acids, characteristic products of ageing of fatty acids with unsaturations in C8 and C9 positions, present in oils.

The C17 fatty acid methyl ester should be indicative of the presence of animal fat. N-alkanes and long chain fatty acids methyl esters are markers of the presence of a wax.